

SUFU

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 53kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SUFU

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

bm

Synonym:

SUFUH; SUFUXL; PRO1280;

Catalog #: A6757

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 51684

Isotype: IgG

Swiss Prot: Q9UMX1

Purity: Affinity purification

For research use only.

Background:

The Hedgehog signaling pathway plays an important role in early human development. The pathway is a signaling cascade that plays a role in pattern formation and cellular proliferation during development. This gene encodes a negative regulator of the hedgehog signaling pathway. Defects in this gene are a cause of medulloblastoma. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).