

TRAPPC10

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 142kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TRAPPC10

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

EHOC1; GT334; TMEM1; TRS30; EHOC-1; TRS130;

Catalog #: A6777

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7109

Isotype: IgG

Swiss Prot: P48553

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a transmembrane protein found in the cis-Golgi complex. The encoded protein is part of the multisubunit transport protein particle (TRAPP) complex and may be involved in vesicular transport from the endoplasmic reticulum to the Golgi. Mutations in this gene could be responsible for the Unverricht-Lundborg type of progressive myoclonus epilepsy, or for autoimmune polyglandular disease type 1.

To place an order, please [Click HERE](#).