

TSC22D3

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 14kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TSC22D3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DIP; GILZ; hDIP; DSIPI; TSC-22R;

Catalog #: A6779

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1831

Isotype: IgG

Swiss Prot: Q99576

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene shares significant sequence identity with the murine TSC-22 and Drosophila shs, both of which are leucine zipper proteins, that function as transcriptional regulators. The expression of this gene is stimulated by glucocorticoids and interleukin 10, and it appears to play a key role in the anti-inflammatory and immunosuppressive effects of this steroid and chemokine. Transcript variants encoding different isoforms have been identified for this gene.

To place an order, please [Click HERE](#).