

UFD1L

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human UFD1L

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

UFD1;

Catalog #: A6783

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7353

Isotype: IgG

Swiss Prot: Q92890

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene forms a complex with two other proteins, nuclear protein localization-4 and valosin-containing protein, and this complex is necessary for the degradation of ubiquitinated proteins. In addition, this complex controls the disassembly of the mitotic spindle and the formation of a closed nuclear envelope after mitosis. Mutations in this gene have been associated with Catch 22 syndrome as well as cardiac and craniofacial defects. Alternative splicing results in multiple transcript variants encoding different isoforms. A related pseudogene has been identified on chromosome 18.

To place an order, please [Click HERE](#).