

XPO5

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 136kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human XPO5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

exp5;

Catalog #: A6790

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57510

Isotype: IgG

Swiss Prot: Q9HAV4

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the karyopherin family that is required for the transport of small RNAs and double-stranded RNA-binding proteins from the nucleus to the cytoplasm. The encoded protein translocates cargo through the nuclear pore complex in a RanGTP-dependent process.

To place an order, please [Click HERE](#).