

ZBED1

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 78kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ZBED1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ALTE; DREF; TRAMP; hDREF;

Catalog #: A6792

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9189

Isotype: IgG

Swiss Prot: O96006

Purity: Affinity purification

For research use only.

Background:

This gene is located in the pseudoautosomal region 1 (PAR1) of X and Y chromosomes. It was earlier identified as a gene with similarity to Ac transposable elements, however, was found not to have transposase activity. Later studies show that this gene product is localized in the nucleus and functions as a transcription factor. It binds to DNA elements found in the promoter regions of several genes related to cell proliferation, such as histone H1, hence may have a role in regulating genes related to cell proliferation. Alternatively spliced transcript variants with different 5' untranslated region have been found for this gene.

To place an order, please [Click HERE](#).