

Bcl9

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 149kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human Bcl9

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

Gm130; 2610202E01Rik; 8030475K17Rik; A330041G23Rik

To place an order, please [Click HERE](#).

Catalog #: A6795

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 607

Isotype: IgG

Swiss Prot: Q9D219

Purity: Affinity purification

For research use only.