

CCS

Reactivity: Human

Tested applications: WB IP

Recommended Dilution: WB 1:500 - 1:2000 IP 1:50 - 1:500

Calculated MW: 29kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human CCS

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A6803

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9973

Isotype: IgG

Swiss Prot: O14618

Purity: Affinity purification

For research use only.

Background:

Copper chaperone for superoxide dismutase specifically delivers Cu to copper/zinc superoxide dismutase and may activate copper/zinc superoxide dismutase through direct insertion of the Cu cofactor.

To place an order, please [Click HERE](#).