

CRY1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200 IF 1:10 - 1:100

Calculated MW: 66kDa

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human CRY1

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PHLL1;

Catalog #: A6807

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1407

Isotype: IgG

Swiss Prot: Q16526

Purity: Affinity purification

For research use only.

Background:

This gene encodes a flavin adenine dinucleotide-binding protein that is a key component of the circadian core oscillator complex, which regulates the circadian clock. This gene is upregulated by CLOCK/ARNTL heterodimers but then represses this upregulation in a feedback loop using PER/CRY heterodimers to interact with CLOCK/ARNTL. Polymorphisms in this gene have been associated with altered sleep patterns. The encoded protein is widely conserved across plants and animals. Loss of the related gene in mouse results in a shortened circadian cycle in complete darkness.

To place an order, please [Click HERE](#).