

DCLRE1B

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 60kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human DCLRE1B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SNM1B; SNMIB; APOLLO

Catalog #: A6808

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 64858

Isotype: IgG

Swiss Prot: Q9H816

Purity: Affinity purification

For research use only.

Background:

DNA interstrand cross-links prevent strand separation, thereby physically blocking transcription, replication, and segregation of DNA. DCLRE1B is one of several evolutionarily conserved genes involved in repair of interstrand cross-links (Dronkert et al., 2000 [PubMed 10848582]).

To place an order, please [Click HERE](#).