

## GTF2H1

---

**Reactivity:** Human Mouse Rat

**Tested applications:** WB

**Recommended Dilution:** WB 1:500 - 1:1000

**Calculated MW:** 62kDa

**Observed MW:** Refer to figures

**Immunogen:**

Recombinant protein of human GTF2H1

**Storage Buffer:**

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

**Synonym:**

P62; BTF2; TFB1; TFIIH;

*To place an order, please [Click HERE](#).*

**Catalog #:** A6814

**Antibody Type:**

Polyclonal Antibody

**Species:** Rabbit

**Gene ID:** 2965

**Isotype:** IgG

**Swiss Prot:** P32780

**Purity:** Affinity purification

For research use only.