

HUS1B

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 31kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human HUS1B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Catalog #: A6817

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 135458

Isotype: IgG

Swiss Prot: Q8NHY5

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is most closely related to HUS1, a component of a cell cycle checkpoint protein complex involved in cell cycle arrest in response to DNA damage. This protein can interact with the check point protein RAD1 but not with RAD9. Overexpression of this protein has been shown to induce cell death, which suggests a related but distinct role of this protein, as compared to the HUS1.

To place an order, please [Click HERE](#).