

NTHL1

Reactivity: Human

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 34kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NTHL1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FAP3; NTH1; OCTS3; hNTH1;

Catalog #: A6820

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 4913

Isotype: IgG

Swiss Prot: P78549

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a DNA N-glycosylase of the endonuclease III family. Like a similar protein in E. coli, the encoded protein has DNA glycosylase activity on DNA substrates containing oxidized pyrimidine residues and has apurinic/apyrimidinic lyase activity.

To place an order, please [Click HERE](#).