

RDM1

Reactivity:Mouse

Tested applications:WB

Recommended Dilution:WB 1:200 - 1:2000

Calculated MW:32kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human RDM1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RAD52B;

Catalog #:A6845

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:201299

Isotype:IgG

Swiss Prot:Q8NG50

Purity:Affinity purification

For research use only.

Background:

This gene encodes a protein involved in the cellular response to cisplatin, a drug commonly used in chemotherapy. The protein encoded by this gene contains two motifs: a motif found in RAD52, a protein that functions in DNA double-strand breaks and homologous recombination, and an RNA recognition motif (RRM) that is not found in RAD52. The RAD52 motif region in RAD52 is important for protein function and may be involved in DNA binding or oligomerization. Alternatively spliced transcript variants encoding different isoforms have been reported.

To place an order, please [Click HERE](#).