

PPP4R3B

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:20 - 1:200

Observed MW: Refer to figures

Immunogen:

A synthetic Peptide of human PPP4R3B

Storage Buffer:

Store at 4. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PSY2; smk1; FLFL2; SMEK2; PP4R3B;

To place an order, please [Click HERE](#).

Catalog #: A6850

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 57223

Isotype: IgG

Swiss Prot: Q5MIZ7

Purity: Affinity purification

For research use only.