

UBE2T

Reactivity: Human Mouse Rat

Tested applications: WB

Recommended Dilution: WB 1:200 - 1:2000

Calculated MW: 23kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human UBE2T

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FANCT; PIG50; HSPC150;

Catalog #: A6853

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 29089

Isotype: IgG

Swiss Prot: Q9NPD8

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene catalyzes the covalent attachment of ubiquitin to protein substrates. Defects in this gene have been associated with Fanconi anemia of complementation group T. Two transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).