

ALOX15

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 75kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ALOX15

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

12-LOX; 15LOX-1; 15-LOX-1;

To place an order, please [Click HERE](#).

Catalog #: A6864

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 246

Isotype: IgG

Swiss Prot: P16050

Purity: Affinity purification

For research use only.