

ALPPL2

Reactivity:Human Mouse

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:57kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human ALPPL2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ALPG; GCAP; ALPPL;

Catalog #:A6866

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:251

Isotype:IgG

Swiss Prot:P10696

Purity:Affinity purification

For research use only.

Background:

There are at least four distinct but related alkaline phosphatases: intestinal, placental, placental-like, and liver/bone/kidney (tissue non-specific). The product of this gene is a membrane bound glycosylated enzyme, localized to testis, thymus and certain germ cell tumors, that is closely related to both the placental and intestinal forms of alkaline phosphatase.

To place an order, please [Click HERE](#).