

C3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 187kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human C3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ASP; C3a; C3b; AHUS5; ARMD9; CPAMD1; HEL-S-62p;

Catalog #: A6879

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 718

Isotype: IgG

Swiss Prot: P01024

Purity: Affinity purification

For research use only.

Background:

Complement component C3 plays a central role in the activation of complement system. Its activation is required for both classical and alternative complement activation pathways. People with C3 deficiency are susceptible to bacterial infection.

To place an order, please [Click HERE](#).