

CKS1B

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 10kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CKS1B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CKS1; ckshs1; PNAS-16; PNAS-18;

Catalog #: A6885

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1163

Isotype: IgG

Swiss Prot: P61024

Purity: Affinity purification

For research use only.

Background:

CKS1B protein binds to the catalytic subunit of the cyclin dependent kinases and is essential for their biological function. The CKS1B mRNA is found to be expressed in different patterns through the cell cycle in HeLa cells, which reflects a specialized role for the encoded protein. At least two transcript variants have been identified for this gene, and it appears that only one of them encodes a protein.

To place an order, please [Click HERE](#).