

DACH1

Reactivity:Human Mouse Rat

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:79kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human DACH1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DACH;

Catalog #:A6895

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:1602

Isotype:IgG

Swiss Prot:Q9UI36

Purity:Affinity purification

For research use only.

Background:

This gene encodes a chromatin-associated protein that associates with other DNA-binding transcription factors to regulate gene expression and cell fate determination during development. The protein contains a Ski domain that is highly conserved from Drosophila to human. Expression of this gene is lost in some forms of metastatic cancer, and is correlated with poor prognosis. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).