

DEFA1

Reactivity: Human Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW: 10kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DEFA1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HP1; MRS; DEF1; HP-1; DEFA2; HNP-1;

Catalog #: A6897

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1667

Isotype: IgG

Swiss Prot: P59665

Purity: Affinity purification

For research use only.

Background:

Defensins are a family of antimicrobial and cytotoxic peptides thought to be involved in host defense. They are abundant in the granules of neutrophils and also found in the epithelia of mucosal surfaces such as those of the intestine, respiratory tract, urinary tract, and vagina. Members of the defensin family are highly similar in protein sequence and distinguished by a conserved cysteine motif. The protein encoded by this gene, defensin, alpha 1, is found in the microbicidal granules of neutrophils and likely plays a role in phagocyte-mediated host defense. Several alpha defensin genes are clustered on chromosome 8. This gene differs from defensin, alpha 3 by only one amino acid. This gene and the gene encoding defensin, alpha 3 are both subject to copy number variation.

To place an order, please [Click HERE](#).