

DLST

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 49kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DLST

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DLTS;

Catalog #: A6901

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 1743

Isotype: IgG

Swiss Prot: P36957

Purity: Affinity purification

For research use only.

Background:

This gene encodes a mitochondrial protein that belongs to the 2-oxoacid dehydrogenase family.

This protein is one of the three components (the E2 component) of the 2-oxoglutarate dehydrogenase complex that catalyzes the overall conversion of 2-oxoglutarate to succinyl-CoA and CO(2). Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).