

FABP6

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 14kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human FABP6

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ILBP; I-15P; I-BAP; ILBP3; ILLBP; I-BABP; I-BALB;

Catalog #: A6906

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2172

Isotype: IgG

Swiss Prot: P51161

Purity: Affinity purification

For research use only.

Background:

This gene encodes the ileal fatty acid binding protein. Fatty acid binding proteins are a family of small, highly conserved, cytoplasmic proteins that bind long-chain fatty acids and other hydrophobic ligands. FABP6 and FABP1 (the liver fatty acid binding protein) are also able to bind bile acids. It is thought that FABPs roles include fatty acid uptake, transport, and metabolism. Transcript variants generated by alternate transcription promoters and/or alternate splicing have been found for this gene.

To place an order, please [Click HERE](#).