

BLOC1S1

Reactivity: Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 17kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human BLOC1S1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RT14; BLOS1; MICoA; GCN5L1;

Catalog #: A6912

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2647

Isotype: IgG

Swiss Prot: P78537

Purity: Affinity purification

For research use only.

Background:

BLOC1S1 is a component of the ubiquitously expressed BLOC1 multisubunit protein complex. BLOC1 is required for normal biogenesis of specialized organelles of the endosomal-lysosomal system, such as melanosomes and platelet dense granules (Starcevic and Dell'Angelica, 2004

To place an order, please [Click HERE](#).