

GIF

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 45kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human GIF

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

IF; INF; IFMH; TCN3;

Catalog #: A6914

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 2694

Isotype: IgG

Swiss Prot: P27352

Purity: Affinity purification

For research use only.

Background:

This gene is a member of the cobalamin transport protein family. It encodes a glycoprotein secreted by parietal cells of the gastric mucosa and is required for adequate absorption of vitamin B12. Vitamin B12 is necessary for erythrocyte maturation and mutations in this gene may lead to congenital pernicious anemia.

To place an order, please [Click HERE](#).