

HBZ

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:16kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human HBZ

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HBZ1; HBZ-T1;

Catalog #:A6920

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:3050

Isotype:IgG

Swiss Prot:P02008

Purity:Affinity purification

For research use only.

Background:

Zeta-globin is an alpha-like hemoglobin. The zeta-globin polypeptide is synthesized in the yolk sac of the early embryo, while alpha-globin is produced throughout fetal and adult life. The zeta-globin gene is a member of the human alpha-globin gene cluster that includes five functional genes and two pseudogenes. The order of genes is: 5' - zeta - pseudozeta - mu - pseudoalpha-1 - alpha-2 -alpha-1 - theta1 - 3'.

To place an order, please [Click HERE](#).