

PPID

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 41kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PPID

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

CYPD; CYP-40;

Catalog #: A6949

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5481

Isotype: IgG

Swiss Prot: Q08752

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the peptidyl-prolyl cis-trans isomerase (PPIase) family. PPIases catalyze the cis-trans isomerization of proline imidic peptide bonds in oligopeptides and accelerate the folding of proteins. This protein has been shown to possess PPIase activity and, similar to other family members, can bind to the immunosuppressant cyclosporin A.

To place an order, please [Click HERE](#).