

RBBP5

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 59kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human RBBP5

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

RBQ3; SWD1;

Catalog #: A6965

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 5929

Isotype: IgG

Swiss Prot: Q15291

Purity: Affinity purification

For research use only.

Background:

This gene encodes a ubiquitously expressed nuclear protein which belongs to a highly conserved subfamily of WD-repeat proteins. The encoded protein binds directly to retinoblastoma protein, which regulates cell proliferation. It interacts preferentially with the underphosphorylated retinoblastoma protein via the E1A-binding pocket B. Three alternatively spliced transcript variants that encode different protein isoforms have been described for this gene.

To place an order, please [Click HERE](#).