

SRSF4

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 57kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SRSF4

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SFRS4; SRP75;

Catalog #: A6978

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6429

Isotype: IgG

Swiss Prot: Q08170

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the arginine/serine-rich splicing factor family. The encoded protein likely functions in mRNA processing.

To place an order, please [Click HERE](#).