

SGCB

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100

Calculated MW: 35kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SGCB

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

A3b; SGC; LGMD2E;

Catalog #: A6979

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6443

Isotype: IgG

Swiss Prot: Q16585

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the sarcoglycan family. Sarcoglycans are transmembrane components in the dystrophin-glycoprotein complex which help stabilize the muscle fiber membranes and link the muscle cytoskeleton to the extracellular matrix. Mutations in this gene have been associated with limb-girdle muscular dystrophy.

To place an order, please [Click HERE](#).