

SSBP1

Reactivity: Human

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 17kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human SSBP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SSBP; mtSSB; Mt-SSB; SOSS-B1;

Catalog #: A6987

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6742

Isotype: IgG

Swiss Prot: Q04837

Purity: Affinity purification

For research use only.

Background:

SSBP1 is a housekeeping gene involved in mitochondrial biogenesis (Tiranti et al., 1995)

To place an order, please [Click HERE](#).