

STK3

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 56kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human STK3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KRS1; MST2;

Catalog #: A6992

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 6788

Isotype: IgG

Swiss Prot: Q13188

Purity: Affinity purification

For research use only.

Background:

This gene encodes a serine/threonine protein kinase activated by proapoptotic molecules indicating the encoded protein functions as a growth suppressor. Cleavage of the protein product by caspase removes the inhibitory C-terminal portion. The N-terminal portion is transported to the nucleus where it homodimerizes to form the active kinase which promotes the condensation of chromatin during apoptosis. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).