

ZNF148

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 89kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ZNF148

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

BERF-1; BFCOL1; ZBP-89; ZFP148; pHZ-52; HT-BETA;

To place an order, please [Click HERE](#).

Catalog #: A7001

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7707

Isotype: IgG

Swiss Prot: Q9UQR1

Purity: Affinity purification

For research use only.