

NELFE

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 43kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NELFE

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

RD; RDP; RDBP; D6S45; NELF-E;

Catalog #: A7006

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 7936

Isotype: IgG

Swiss Prot: P18615

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is part of a complex termed negative elongation factor (NELF) which represses RNA polymerase II transcript elongation. This protein bears similarity to nuclear RNA-binding proteins; however, it has not been demonstrated that this protein binds RNA. The protein contains a tract of alternating basic and acidic residues, largely arginine (R) and aspartic acid (D). The gene localizes to the major histocompatibility complex (MHC) class III region on chromosome 6.

To place an order, please [Click HERE](#).