

NPFF

Reactivity:Human

Tested applications:WB IHC

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW:12kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human NPFF

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

FMRFAL;

Catalog #:A7020

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:8620

Isotype:IgG

Swiss Prot:O15130

Purity:Affinity purification

For research use only.

Background:

FMFRamide-related protein precursor plays a role in the regulation of heart rate and blood pressure and the modulation of morphine-induced antinociception. FMRFAL encodes a preproprotein which is cleaved to form two active peptides with similar function.

To place an order, please [Click HERE](#).