

ASMTL

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 69kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ASMTL

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

f

Synonym:

ASTML; ASMTLX; ASMTLY;

Catalog #: A7021

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 8623

Isotype: IgG

Swiss Prot: O95671

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene has an N-terminus that is similar to the multicopy associated filamentation (maf) protein of Bacillus subtilis and to orfE of Escherichia coli, while the C-terminus is similar to N-acetylserotonin O-methyltransferase. This gene is located in the pseudoautosomal region 1 (PAR1) of X and Y chromosomes. Three transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).