

EFTUD2

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:200

Calculated MW: 109kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human EFTUD2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Concentration:

j

Synonym:

MFDM; MFDGA; Snu114; Snrp116; SNRNP116; U5-116KD;

Catalog #: A7040

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 9343

Isotype: IgG

Swiss Prot: Q15029

Purity: Affinity purification

For research use only.

Background:

This gene encodes a GTPase which is a component of the spliceosome complex which processes precursor mRNAs to produce mature mRNAs. Mutations in this gene are associated with mandibulofacial dysostosis with microcephaly. Multiple transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).