

DPP3

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 83kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human DPP3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DPPIII;

Catalog #: A7054

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10072

Isotype: IgG

Swiss Prot: Q9NY33

Purity: Affinity purification

For research use only.

Background:

This gene encodes a protein that is a member of the M49 family of metallopeptidases. This cytoplasmic protein binds a single zinc ion with its zinc-binding motif (HELLGH) and has post-proline dipeptidyl aminopeptidase activity, cleaving Xaa-Pro dipeptides from the N-termini of proteins. Increased activity of this protein is associated with endometrial and ovarian cancers. Alternatively spliced transcript variants have been found for this gene.

To place an order, please [Click HERE](#).