

PIAS3

Reactivity: Human Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:100 IF 1:50 - 1:100

Calculated MW: 68kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PIAS3

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

ZMIZ5;

Catalog #: A7060

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10401

Isotype: IgG

Swiss Prot: Q9Y6X2

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the PIAS

To place an order, please [Click HERE](#).