

VTI1B

Reactivity:Human Mouse

Tested applications:WB IHC IF

Recommended Dilution:WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW:27kDa

Observed MW:Refer to figures

Immunogen:

Recombinant protein of human VTI1B

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

VTI1; VTI2; VTI1L; v-SNARE; vti1-rp1; VTI1-LIKE;

To place an order, please [Click HERE](#).

Catalog #:A7062

Antibody Type:

Polyclonal Antibody

Species:Rabbit

Gene ID:10490

Isotype:IgG

Swiss Prot:Q9UEU0

Purity:Affinity purification

For research use only.