

ENOX2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 70kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ENOX2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

APK1; tNOX; COVA1;

Catalog #: A7063

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10495

Isotype: IgG

Swiss Prot: Q16206

Purity: Affinity purification

For research use only.

Background:

This gene is a tumor-specific member of the ECTO-NOX family of genes that encode cell surface NADH oxidases. The encoded protein has two enzymatic activities: catalysis of hydroquinone or NADH oxidation, and protein disulfide interchange. The protein also displays prion-like properties. Alternative splicing results in multiple transcript variants.

To place an order, please [Click HERE](#).