

PLK2

Reactivity: Human Mouse

Tested applications: WB

Recommended Dilution: WB 1:500 - 1:2000

Calculated MW: 78kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PLK2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

SNK; hSNK; hPlk2;

Catalog #: A7066

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 10769

Isotype: IgG

Swiss Prot: Q9NYY3

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the polo family of serine/threonine protein kinases that have a role in normal cell division. This gene is most abundantly expressed in testis, spleen and fetal tissues, and its expression is inducible by serum, suggesting that it may also play an important role in cells undergoing rapid cell division. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.

To place an order, please [Click HERE](#).