

TUSC2

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 12kDa

Observed MW: Refer to Figures

Immunogen:

Recombinant protein of human TUSC2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

PAP; FUS1; PDAP2; C3orf11;

Catalog #: A7075

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 11334

Isotype: IgG

Swiss Prot: O75896

Purity: Affinity purification

For research use only.

Background:

This gene is a highly conserved lung cancer candidate gene. No other information about this gene is currently available.

To place an order, please [Click HERE](#).