

TRIM32

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 72kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human TRIM32

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

HT2A; BBS11; TATIP; LGMD2H;

Catalog #: A7079

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 22954

Isotype: IgG

Swiss Prot: Q13049

Purity: Affinity purification

For research use only.

Background:

The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to cytoplasmic bodies. The protein has also been localized to the nucleus, where it interacts with the activation domain of the HIV-1 Tat protein. The Tat protein activates transcription of HIV-1 genes.

To place an order, please [Click HERE](#).