

ICOSLG

Reactivity: Human Mouse Rat

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 33kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human ICOSLG

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

B7H2; GL50; B7-H2; B7RP1; CD275; ICOSL; LICOS; B7RP-1; ICOS-L;

To place an order, please [Click HERE](#).

Catalog #: A7080

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23308

Isotype: IgG

Swiss Prot: O75144

Purity: Affinity purification

For research use only.