

CLASP1

Reactivity: Human Mouse Rat

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 169kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human CLASP1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

MAST1;

Catalog #: A7081

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23332

Isotype: IgG

Swiss Prot: Q7Z460

Purity: Affinity purification

For research use only.

Background:

CLASPs, such as CLASP1, are nonmotor microtubule-associated proteins that interact with CLIPs (e.g., CLIP170; MIM 179838). CLASP1 is involved in the regulation of microtubule dynamics at the kinetochore and throughout the spindle (Maiato et al., 2003)

To place an order, please [Click HERE](#).