

PRND

Reactivity: Human

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 20kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human PRND

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

DPL; PrPLP; DOPPEL; dJ1068H6.4;

Catalog #: A7083

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 23627

Isotype: IgG

Swiss Prot: Q9UKY0

Purity: Affinity purification

For research use only.

Background:

This gene is found on chromosome 20, approximately 20 kbp downstream of the gene encoding cellular prion protein, to which it is biochemically and structurally similar. The protein encoded by this gene is a membrane glycosylphosphatidylinositol-anchored glycoprotein that is found predominantly in testis. Mutations in this gene may lead to neurological disorders.

To place an order, please [Click HERE](#).