

NTMT1

Reactivity: Human Mouse

Tested applications: WB IHC IF

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200 IF 1:50 - 1:100

Calculated MW: 25kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human NTMT1

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

NRMT; NTM1A; AD-003; HONT1A; C9orf32; METTL11A;

Catalog #: A7098

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 28989

Isotype: IgG

Swiss Prot: Q9BV86

Purity: Affinity purification

For research use only.

Background:

The METTL11A gene encodes an N-terminal methyltransferase for the RAN (MIM 601179) guanine nucleotide exchange factor regulator of chromosome condensation 1 (RCC1; MIM 179710). METTL11A enzyme alpha-N-methylates other protein targets such as SET (MIM 600960) and RB (MIM 180200).

To place an order, please [Click HERE](#).