

KCNIP2

Reactivity: Human Mouse

Tested applications: WB IHC

Recommended Dilution: WB 1:500 - 1:2000 IHC 1:50 - 1:200

Calculated MW: 31kDa

Observed MW: Refer to figures

Immunogen:

Recombinant protein of human KCNIP2

Storage Buffer:

Store at -20. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Synonym:

KCHIP2;

Catalog #: A7100

Antibody Type:

Polyclonal Antibody

Species: Rabbit

Gene ID: 30819

Isotype: IgG

Swiss Prot: Q9NS61

Purity: Affinity purification

For research use only.

Background:

This gene encodes a member of the family of voltage-gated potassium (Kv) channel-interacting proteins (KCNIPs), which belongs to the recoverin branch of the EF-hand superfamily. Members of the KCNIP family are small calcium binding proteins. They all have EF-hand-like domains, and differ from each other in the N-terminus. They are integral subunit components of native Kv4 channel complexes. They may regulate A-type currents, and hence neuronal excitability, in response to changes in intracellular calcium. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified from this gene.

To place an order, please [Click HERE](#).